


# iCup<sup>A.D.</sup> Drug Screen<sup>®</sup>



## iCup<sup>®</sup> Drug Screen

*Corporate Drug Testing Made Simple.*

- Room temperature storage
- Built-in validity test
- Detects 3-10 drugs
- 24-month shelf life
- Simple procedure
- Fast results

**Brands You Can Trust<sup>™</sup>**




# iCup<sup>®</sup> A.D. Drug Screen<sup>®</sup>

## iCup<sup>®</sup> Drug Screen Corporate Drug Testing Made Simple.

Detect the presence of drugs and adulterants in minutes. It's drug testing efficiency you can count on.

### Integrated testing solution:

- iCup<sup>®</sup> with temperature strip
- Security Seal
- Two-part result form available
- Photocopy Template




iCup<sup>®</sup> AD 5-panel shown as a negative result.

Drug Name:	Abbreviation:	Cutoff Levels:
Amphetamine.....	AMP.....	1000ng/ml
Barbiturates.....	BAR.....	300ng/ml
Benzodiazepines.....	BZO.....	300ng/ml
Cocaine.....	COC.....	300ng/ml
Ecstasy.....	MDMA.....	500ng/ml
Marijuana.....	THC.....	50ng/ml
Methadone.....	MTD.....	300ng/ml
Methamphetamine.....	mAMP.....	1000ng/ml
Morphine.....	MOR.....	300ng/ml
Opiates.....	OPI.....	2000ng/ml
Oxycodone.....	OXY.....	100ng/ml
Phencyclidine.....	PCP.....	25ng/ml
Propoxyphene.....	PPX.....	300ng/ml
Tricyclic Antidepressants.....	TCA.....	1000ng/ml

## Instructions


### 1 Collect Specimen




### 2 Peel Back Label



### 3 Read the Results



Read results at five (5) minutes.  
Results stable up to sixty (60) minutes.



SAMPLE ONLY - Color chart NOT Exact